

E 130
Insert
for the optimum loading of 10 trays.



- equipped with 11 holders, height 170 mm, distance apart 35 mm
- For use in lower basket

EAN: 4002511129375 / Material number: 02648950 / Old Material Number: 69513001D

Product category	
Thermal disinfectors	•
Large-chamber thermal disinfectors, dental	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for trays	•
Area of application	
Reprocessing of trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	35
Rail height in mm	170
Capacity	
Trays [quantity]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	180
External dimensions, net width in mm	180
External dimensions, net depth in mm	445
External dimensions, gross height in mm	210
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	470
Net weight in kg	0.76
Gross weight in kg	0.90