


O 176

Upper basket

for the optimum loading of instruments for MIS, arthroscopy and urology.



- Equipped with 14 injector nozzles/funnels
- Left inserts, right hollow instruments (up to 500 mm)
- Vertical clearance left is 360 +/- 20 mm
- Powder dispensing not possible
- No spray arm, with connection for hot air drying

EAN: 4002511114890 / Material number: 03807420 / Old Material Number: 69117601D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
<b>Capacity</b>	
Lumened instruments [number]	10
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	409
External dimensions, net width in mm	475
External dimensions, net depth in mm	531
External dimensions, gross height in mm	545
External dimensions, gross width in mm	400
External dimensions, gross depth in mm	530
Net weight in kg	4.80
Gross weight in kg	6.50
<b>Companion pack</b>	
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	2
E 448 silicone hose with female Luer Lock adapter [number]	2
E 452 injector nozzle 2.5 x 60 mm [number]	1
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	1
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
Spray jet, Ø 2.5 mm [number]	2
Spray jet, Ø 4 mm [number]	4
Funnel [number]	4