


O 183

Upper basket

for the optimum loading of instruments for MIS, arthroscopy and urology.



- Left side: 10 injector nozzles or free accommodation of inserts
- Right side: 14 injector nozzles/ funnels for hollow instruments
- Vertical clearance 285 +/- 30mm; hollow instruments 370 +/- 30mm
- Powder dispensing not possible

EAN: 4002511114920 / Material number: 03807470 / Old Material Number: 69118301D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
<b>Capacity</b>	
Lumened instruments [number]	14
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	270
External dimensions, net width in mm	475
External dimensions, net depth in mm	530
External dimensions, gross height in mm	410
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	555
Net weight in kg	4.00
Gross weight in kg	4.80