


O 184

Upper basket

for optimum loading of centrifuge tubes, phials or test tubes, for example.



- Equipped with 96 injector nozzles
- Equipped with connection for hot air drying
- Vertical clearance: 170 mm
- Powder dispensing not possible

EAN: 4002511114937 / Material number: 03807480 / Old Material Number: 69118401D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Basket for vials	•
Basket for test tubes	•
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	170
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•
Capacity	
Test tubes [number]	96
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	468
External dimensions, net width in mm	475
External dimensions, net depth in mm	531
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	4.50
Gross weight in kg	5.70
Companion pack	
ID 110, injector nozzle, Ø 2.5 mm, length 110 mm [number]	96