

U 167

Lower basket

for the optimum loading of anaesthetic utensils.



- Suitable for 4 soda lime containers, 9 secretion flasks & various utensils
- 20 holders, height 200 mm mm, spacing approx. 95 mm mm

EAN: 4002511114784 / Material number: 03807600 / Old Material Number: 69116702D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Lower basket	•
Mobile unit for anaesthetic instruments	•
<b>Product properties</b>	
Material	Wire, rilsan-coated
Colour	White
Distance between rails in mm	95
Rail height in mm	200
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	220
External dimensions, net width in mm	535
External dimensions, net depth in mm	516
External dimensions, gross height in mm	240
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	3.10
Gross weight in kg	4.60