


E 111

Insert

For the optimum loading of beakers from 600 to 1,000 ml.



- Suitable for 8 glass beakers between 600 – 1000 ml
- equipped with 8 x 3 stay bars
- For use in lower basket

EAN: 4002511115743 / Material number: 03808420 / Old Material Number: 69511101D

| Product classification | |
|--|--------------------------------------|
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| G 7825 | • |
| PG 8527 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| PLW 6011 | • |
| PLW 6111 | • |
| PLW 8505 | • |
| Product category | |
| Insert/module for beakers | • |
| Product properties | |
| Material | Stainless steel, Wire, rilsan-coated |
| Colour | Stainless Steel, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Beakers [quantity] | 8 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 205 |
| External dimensions, net width in mm | 230 |
| External dimensions, net depth in mm | 460 |
| External dimensions, gross height in mm | 230 |
| External dimensions, gross width in mm | 250 |
| External dimensions, gross depth in mm | 480 |
| Net weight in kg | 0.96 |
| Gross weight in kg | 1.10 |