

E 124

Insert

For the optimum loading of infusion flasks. Surface area 100 x 100 mm.



- 16 compartments, dimension: 100 x 100 mm (bottle), 48 x 48 mm (neck)
- For use in lower basket

EAN: 4002511115811 / Material number: 03808500 / Old Material Number: 69512401D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	16
Area of application	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	148
External dimensions, net width in mm	460
External dimensions, net depth in mm	460
External dimensions, gross height in mm	170
External dimensions, gross width in mm	480
External dimensions, gross depth in mm	480
Net weight in kg	1.44
Gross weight in kg	1.80