


E 125

Insert

For the optimum loading of infusion flasks. Surface area 125 x 125 mm.



- 9 compartments, dimension: 125 x 125 mm (bottle), 55 x 55 mm (neck)
- For use in lower basket

EAN: 4002511115828 / Material number: 03808510 / Old Material Number: 69512501D

Product classification	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	9
Area of application	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	224
External dimensions, net width in mm	460
External dimensions, net depth in mm	460
External dimensions, gross height in mm	245
External dimensions, gross width in mm	480
External dimensions, gross depth in mm	480
Net weight in kg	1.81
Gross weight in kg	2.56