

E 128

Insert

For the optimum loading of infusion flasks. Surface area 71 x 71 mm.



- 24 compartments, dimension: 71 x 71 mm (bottle), 46 x 46 mm (neck)

EAN: 4002511115859 / Material number: 03808540 / Old Material Number: 69512801D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	24
Area of application	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	103
External dimensions, net width in mm	445
External dimensions, net depth in mm	445
External dimensions, gross height in mm	160
External dimensions, gross width in mm	475
External dimensions, gross depth in mm	470
Net weight in kg	1.49
Gross weight in kg	1.95