

E 135/1

Insert

for optimum loading of 19 baby bottles up to 110 ml.



- Suitable for 19 baby bottles x 110 ml
- Suitable for bottles sized 51 x 51 mm, neck: 45 x 45 mm
- For upper basket A 102 and lower basket A 151

EAN: 4002511115910 / Material number: 03808640 / Old Material Number: 69513502D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of baby bottles	•
Capacity	
Baby bottles [number]	19
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	135
External dimensions, net width in mm	192
External dimensions, net depth in mm	447
External dimensions, gross height in mm	195
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	465
Net weight in kg	2.02
Gross weight in kg	2.10