

E 351

Injector nozzle

For the optimum cleaning of the wash items. Size 4 x 160 mm.



- Suitable for use in a mobile injector unit
- Can be combined with clip E 353
- Ø 4 x 160 mm, screw connection

EAN: 4002511117075 / Material number: 03809500 / Old Material Number: 69735101D

| | |
|--|-----------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| G 7825 | • |
| PG 8527 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| Product category | |
| Injector nozzle | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 160 |
| External dimensions, net width in mm | 4 |
| External dimensions, net depth in mm | 4 |
| External dimensions, gross height in mm | 12 |
| External dimensions, gross width in mm | 80 |
| External dimensions, gross depth in mm | 120 |
| Diameter in mm | 4 |
| Net weight in kg | 0.01 |
| Gross weight in kg | 0.02 |