


ID 140
 Injector nozzle
 with plastic supports, Ø 4 mm, length 140 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring flasks and test tubes, for example

EAN: 4002511116986 / Material number: 03810340 / Old Material Number: 69714001D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	55
External dimensions, net depth in mm	140
External dimensions, gross height in mm	55
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	145
Diameter in mm	4
Net weight in kg	0.01
Gross weight in kg	0.02