

E 137

Insert

for optimum loading of 56 half petri dishes.



- Suitable for 56 half petri dishes with Ø 100 mm
- equipped with 56 holders, height 70 mm, distance apart approx. 26 mm
- For use in lower basket
- For use with insert E 136

EAN: 4002511115958 / Material number: 03830290 / Old Material Number: 69513701D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert/module for half petri dishes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	26
Rail height in mm	70
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Half petri dishes [number]	56
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	485
External dimensions, net width in mm	445
External dimensions, net depth in mm	95
External dimensions, gross height in mm	140
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	580
Net weight in kg	2.02
Gross weight in kg	3.45