

E 142

Insert

For holding different instruments.



- 1 mm wire, 5 mm mesh size, 5 mm frames
- equipped with 2 hinged handles
- Maximum load 10 kg

EAN: 4002511115972 / Material number: 03830300 / Old Material Number: 69514201D

| Product category                            |                 |
|---|-----------------|
| Washer-disinfectors, laboratory             | •               |
| Washer-disinfectors, medical                | •               |
| Large-capacity washer-disinfectors, medical | •               |
| Machine affiliation                         |                 |
| G 7824                                      | •               |
| G 7826                                      | •               |
| PG 8527                                     | •               |
| PG 8528                                     | •               |
| G 7831                                      | •               |
| PG 8535                                     | •               |
| PG 8536                                     | •               |
| PG 8581                                     | •               |
| PG 8562                                     | •               |
| PG 8582                                     | •               |
| PG 8582 CD                                  | •               |
| PG 8583                                     | •               |
| PG 8583 CD                                  | •               |
| PG 8591                                     | •               |
| PG 8592                                     | •               |
| PG 8593                                     | •               |
| Product category                            |                 |
| Mesh tray                                   | •               |
| Product properties                          |                 |
| Material                                    | Stainless steel |
| Colour                                      | Stainless Steel |
| Area of application                         |                 |
| Reprocessing of laboratory glassware        | •               |
| Reprocessing of GYN instruments             | •               |
| Reprocessing of ENT instruments             | •               |
| Reprocessing of surgical instruments        | •               |
| Standard electrical connection              |                 |
| Number of phases                            | 0               |
| Voltage in V                                | 0               |
| Frequency in Hz                             | 0               |
| Total rated load in kW                      | 0.00            |
| Fuse rating in A                            | 0               |
| Dimensions and weight                       |                 |
| External dimensions, net height in mm       | 55              |
| External dimensions, net width in mm        | 255             |
| External dimensions, net depth in mm        | 480             |
| External dimensions, gross height in mm     | 75              |
| External dimensions, gross width in mm      | 280             |
| External dimensions, gross depth in mm      | 515             |
| Net weight in kg                            | 0.70            |
| Gross weight in kg                          | 0.86            |