

MT  
Trolleys  
for loading & unloading washer-&disinfectors.



- For the storage and transport of baskets and inserts
- Divided into 4 levels that can be adjusted by 102.5 mm
- Moveable with 4 lockable castors
- Loading volume 549 x 599 mm

EAN: 4002511820005 / Material number: 03893140 / Old Material Number: 69200102D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Storage/transport system	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	1985
External dimensions, net width in mm	662
External dimensions, net depth in mm	616
External dimensions, gross height in mm	2030
External dimensions, gross width in mm	680
External dimensions, gross depth in mm	680
Net weight in kg	23.00
Gross weight in kg	29.00