


A 14

Lid

For fixing delicate wash items.



- Lid in stainless steel
- 7 x 7 mm perforations, 3 mm slat

EAN: 4002511860995 / Material number: 03981970 / Old Material Number: 69745002D

Product classification	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Lid	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	20
External dimensions, net width in mm	210
External dimensions, net depth in mm	210
External dimensions, gross height in mm	30
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	250
Net weight in kg	0.20
Gross weight in kg	0.29