


A 16

Lid

For fixing delicate wash items.



- 1 mm wire, 5 mm mesh size, 4 mm frames

EAN: 4002511860971 / Material number: 03981990 / Old Material Number: 69745004D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Lid	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	19
External dimensions, net width in mm	230
External dimensions, net depth in mm	215
External dimensions, gross height in mm	40
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	255
Net weight in kg	0.15
Gross weight in kg	0.30