

E 451

Mesh trays

For optimum accommodation of different instruments.



- Mesh tray with lid for small pieces
- Wire mesh: Base 1 mm, sides 0.8 mm, lid 1 mm
- Mesh width: Base 3 mm, sides 1.7 mm, lid 3 mm

EAN: 4002511908154 / Material number: 04203910 / Old Material Number: 69545101D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Mesh tray	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	55
External dimensions, net width in mm	160
External dimensions, net depth in mm	240
External dimensions, gross height in mm	75
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	245
Net weight in kg	0.79
Gross weight in kg	0.85