

E 451

Mesh trays

For optimum accommodation of different instruments.



- Mesh tray with lid for small pieces
- Wire mesh: Base 1 mm, sides 0.8 mm, lid 1 mm
- Mesh width: Base 3 mm, sides 1.7 mm, lid 3 mm

EAN: 4002511908154 / Material number: 04203910 / Old Material Number: 69545101D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	160
External dimensions, net depth in mm	240
External dimensions, gross height in mm	75
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	245
Net weight in kg	0.79
Gross weight in kg	0.85