

E 454

Injector nozzle

for optimum cleaning of different wash items, such as trocar sleeves.



- Suitable for trocar sleeves with Ø 10-15 mm, 8 x 150 mm
- Equipped with a spring for height adjustment

EAN: 4002511908512 / Material number: 04203940 / Old Material Number: 69745401D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net depth in mm	150
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Diameter in mm	8
Net weight in kg	0.03
Gross weight in kg	0.04