

E 458

Insert

for the optimum loading of screw-cap teats.



- Suitable for 36 screw connection teats
- 36 compartments at 60 x 60
- equipped with a flip-up and lockable lid

EAN: 4002511919358 / Material number: 04235200 / Old Material Number: 69545801D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	36
<b>Area of application</b>	
Reprocessing of baby bottles	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	63
External dimensions, net width in mm	215
External dimensions, net depth in mm	445
External dimensions, gross height in mm	450
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	220
Net weight in kg	1.50
Gross weight in kg	1.60