

E 147/1

Insert

for the optimum loading of 10-12 mouthwash glasses.



- For 10-12 mouth rinse tumblers, max. Ø 80 mm
- Plastic-coated

EAN: 4002512122313 / Material number: 04746620 / Old Material Number: 68514703D

<b>Product category</b>	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Large-chamber thermal disinfectors, dental	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Other inserts	•
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Product properties</b>	
Material	Wire, rilsan-coated
Colour	White
<b>Capacity</b>	
Mouth rinsing cups [number]	12
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	155
External dimensions, net width in mm	220
External dimensions, net depth in mm	455
External dimensions, gross height in mm	260
External dimensions, gross width in mm	510
External dimensions, gross depth in mm	240
Net weight in kg	1.29
Gross weight in kg	1.59