

E 466

Injector nozzle

for the optimum loading of breathing bags 8 x 333 mm.



EAN: 4002512223744 / Material number: 04966620 / Old Material Number: 69746601D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net depth in mm	333
External dimensions, gross height in mm	140
External dimensions, gross width in mm	260
External dimensions, gross depth in mm	150
Diameter in mm	8
Net weight in kg	0.85
Gross weight in kg	1.00