

E 485  
Bracket  
for the optimum loading of 3 kidney trays.



- Several holders will fit in the E 484 insert

EAN: 4002512398725 / Material number: 05109210 / Old Material Number: 69548501D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Insert for kidney dishes/basins	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of kidney dishes	•
<b>Capacity</b>	
Kidney dishes [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	195
External dimensions, net width in mm	80
External dimensions, net depth in mm	464
External dimensions, gross height in mm	70
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	505
Net weight in kg	0.41
Gross weight in kg	0.56