


E 488  
Bracket  
for the optimum loading of respiratory masks.



- Suitable for 3 respirator masks

EAN: 4002512398756 / Material number: 05109260 / Old Material Number: 69548801D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of different utensils	•
<b>Capacity</b>	
Breathing masks [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	31
External dimensions, net width in mm	240
External dimensions, net depth in mm	210
External dimensions, gross height in mm	75
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	510
Net weight in kg	0.27
Gross weight in kg	0.39