

DOS 10/30
Dispenser pump
For dispensing liquid media.



- For liquid media, 10 ml/30 seconds
- Complete with hoses and siphon (330 mm) for 10 l canisters

EAN: 4002512417549 / Material number: 05267410 / Old Material Number: 69225002D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7825	•
G 7826	•
Product category	
Dispenser system	•
Product properties	
Length of suction lance in mm	333
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	116
External dimensions, net width in mm	63
External dimensions, net depth in mm	102
External dimensions, gross height in mm	300
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	425
Net weight in kg	4.50
Gross weight in kg	5.00