

DOS 60/30

Dispenser pump

For dispensing liquid neutral/alkaline media.



- For dispensing liquid cleaning agent
- For dispensing chemical disinfectant
- For retro-fitting by Miele Service

EAN: 4002512417532 / Material number: 05267420 / Old Material Number: 69225003D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7825	•
G 7826	•
Product category	
Dispenser system	•
Product properties	
Production rate in ml/60 s	120
Length of suction lance in mm	333
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	310
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	430
Net weight in kg	1.70
Gross weight in kg	1.95