

E 476

Holder

for micro instruments with a Ø of 4-8 mm.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002513239454 / Material number: 05692400 / Old Material Number: 69747601D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Holder	•
<b>Product properties</b>	
Material	Plastic
Colour	Grey
<b>Area of application</b>	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.07
Gross weight in kg	0.08