

DOS G 60

Dispenser module for 5-10 litre canisters

With level indicator for the automatic dispensing of liquid cleaning agents.



- For use with PG 8535
- Adjustable via electronic machine controls
- 300 mm siphon for 5 and 10 litre containers
- Length of connection cable: 1.9 m

EAN: 4002513252255 / Material number: 05703990 / Old Material Number: 69747005D

Product classification	
Fresh-water dishwashers	•
Tank dishwashers	•
Washer-disinfectors, laboratory	•
Machine affiliation	
G 8066	•
Product category	
External dispensing system	•
Product properties	
Material	Plastic
Colour	Grey
Production rate in ml/60 s	120
Length of suction lance in mm	333
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	230-240
Frequency in Hz	50
Total rated load in kW	0.02
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	117
External dimensions, net width in mm	235
External dimensions, net depth in mm	107
External dimensions, gross height in mm	295
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	420
Net weight in kg	2.60
Gross weight in kg	2.78