

E 969
Insert
for holding different utensils.



- Perforated plate 7 x 7 x 3 mm
- Suitable for E 900-4/1, E 935/1, E 975/2, E 941 with E 945*

* In E 941 only suitable with E 945.

EAN: 4002513298765 / Material number: 05746240 / Old Material Number: 69596901D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
Product category	
Insert/module for beakers	•
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	122
External dimensions, net width in mm	363
External dimensions, net depth in mm	533
External dimensions, gross height in mm	165
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	535
Net weight in kg	1.33
Gross weight in kg	2.23