


E 965
Insert
for the optimum loading of beakers 250 - 600 ml.



- Equipped with 15 x 3 stay bars for beakers from 250 ml to 600 ml

EAN: 4002513369229 / Material number: 05848290 /
Old Material Number: 69596501D

| | |
|--|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PG 8527 | • |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Insert/module for beakers | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Beakers [quantity] | 98 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 173 |
| External dimensions, net width in mm | 357 |
| External dimensions, net depth in mm | 522 |
| External dimensions, gross height in mm | 280 |
| External dimensions, gross width in mm | 540 |
| External dimensions, gross depth in mm | 560 |
| Net weight in kg | 2.07 |
| Gross weight in kg | 2.86 |