

E 960/1

Insert

for the optimum loading of wide necked Erlenmeyer flasks.



- Equipped with 20 large and 26 small spring hooks
- Suitable for Erlenmeyer flasks & measuring cylinders

EAN: 4002513400953 / Material number: 05892360 / Old Material Number: 69596002D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
<b>Product category</b>	
Insert/module for wide-necked glassware	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Erlenmeyer flasks [number]	46
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	185
External dimensions, net width in mm	357
External dimensions, net depth in mm	522
External dimensions, gross height in mm	280
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	560
Net weight in kg	0.80
Gross weight in kg	0.90