


E 450/1
Mobile unit
for the optimum loading of MIS instruments.




- equipped with 42 holding options on 2 levels
- Equipped with connection for hot air drying
- Loading options: E 451, E 457, E 460, E 473, E 444
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002513505504 / Material number: 06059490 / Old Material Number: 69545003D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit for MIS instruments	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of MIS instruments	•
Capacity	
Lumened instruments [number]	42
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	515
External dimensions, gross height in mm	500
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	530
Net weight in kg	9.50
Gross weight in kg	11.06
Companion pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	15
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	5
E 445 – 12 covers for E 442 irrigation sleeve, opening with Ø 6 mm [number]	1
E 446 – 12 covers for E 443 irrigation sleeve, opening with Ø 10 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	6
E 449 male adapter, for female Luer Lock [number]	5
E 451 insert 1/6 mesh tray [number]	1
E 452 injector nozzle 2.5 x 60 mm [number]	3
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	8
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	6
E 456 spring stay for MIS instruments [number]	4

E 450/1
Mobile unit
for the optimum loading of MIS instruments.

EAN: 4002513505504 / Material number: 06059490 / Old Material Number: 69545003D



E 464 holder for injector nozzle E 454 [number]	3
E 472 spring clip for injector nozzles Ø 4 mm [number]	2