

E 427

Insert

for the optimum loading of 6 laryngoscopes.



- Suitable for 6 laryngoscopes

EAN: 4002513808209 / Material number: 06442120 / Old Material Number: 69542701D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert/module for anaesthetic instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	92
External dimensions, net width in mm	210
External dimensions, net depth in mm	134
External dimensions, gross height in mm	164
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	250
Net weight in kg	0.90
Gross weight in kg	1.00