

E 461/2

Mobile unit

for the optimum cleaning of breathing tubes.



- Suitable for 12 breathing tubes up to 1.5m long, spiral storage
- Suitable for: E 432, E 433, E 434, E 430/1, A3
- Equipped with connection for hot air drying
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002513984057 / Material number: 06641180 / Old Material Number: 69546103D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit for anaesthetic instruments	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of anaesthetic instruments	•
Capacity	
Hoses [number]	12
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	515
External dimensions, gross height in mm	546
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	515
Net weight in kg	8.77
Gross weight in kg	10.80
Companion pack	
E 430/1 1/3 mesh tray [number]	1
E 432 holder for concertina breathing tubes [number]	1
E 433 holder for silicone breathing tubes [number]	3
E 434 holder for paediatric breathing tubes [number]	1