

E 435/3

Mobile unit

for optimum loading of approx. 2 anaesthesia sets.



- Suitable for about 2 anaesthesia sets.
- 6 nozzles with spring holders for breathing tubes up to 1.5 m
- Storage on spiral
- Equipped with connection for hot air drying
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002513984064 / Material number: 06641190 / Old Material Number: 69543504D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
Product category	
Mobile unit for anaesthetic instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Capacity	
AN sets [number]	2
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	507
External dimensions, net width in mm	535
External dimensions, net depth in mm	515
External dimensions, gross height in mm	510
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	520
Net weight in kg	6.70
Gross weight in kg	7.00
Companion pack	
E 430/1 1/3 mesh tray [number]	1
E 432 holder for concertina breathing tubes [number]	1