


E 103/1

Insert

For the optimum loading of test tubes up to 12 x 75 mm.



- Suitable for about 200 Test tubes (up to 12 x 75 mm)
- Divided into 6 surfaces
- The mesh width is 8 x 8 mm
- Including lid A 13

EAN: 4002514231051 / Material number: 06907630 / Old Material Number: 69510302D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	200
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	122
External dimensions, net width in mm	200
External dimensions, net depth in mm	220
External dimensions, gross height in mm	140
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	250
Net weight in kg	0.54
Gross weight in kg	0.74
<b>Companion pack</b>	
A 13 lid [number]	1