

E 907/1

Mesh trays

for the optimum loading of different small-scale instruments.



- With lid for small parts
- Mesh width 3 x 1 mm
- Hook for hanging in E 905/1 and E 903/1

EAN: 4002514256320 / Material number: 06943820 / Old Material Number: 69590702D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	48
External dimensions, net width in mm	129
External dimensions, net depth in mm	170
External dimensions, gross height in mm	70
External dimensions, gross width in mm	155
External dimensions, gross depth in mm	230
Net weight in kg	0.30
Gross weight in kg	0.37