

TK 180
Tumble dryer basket
For hanging in the dryer.



- Suitable for PT 5186 and PT 7186
- For drying delicate textiles
- Ideal for drying sports shoes for example

EAN: 4002514715698 / Material number: 07289550 / Old Material Number: 59018001D

Product classification	
Tumble dryers	•
Machine/accessory category	
PT 5186	•
PT 7186	•
PT 7188	•
PWT 6089	•
Product category	
Dryer basket	•
Product properties	
Electrical connection	OHNE
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	475
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight [kg]	0.96
Gross weight in kg	2.94