

O 191/2

Upper basket

for the optimum loading of mesh trays.



- Vertical clearance: 115 +/- 20 mm
- Loading width: 480 mm
- Loading depth: 450 mm
- With integrated magnetic spray arm above the basket
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556192 / Material number: 07333060 / Old Material Number: 69119102D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Carrier/basket carrier	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
DIN mesh trays [number]	2
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	180
External dimensions, net width in mm	531
External dimensions, net depth in mm	475
External dimensions, gross height in mm	540
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	525
Net weight in kg	2.40
Gross weight in kg	3.30