

ASA 1

Magnetic spray arm

for spray arm monitoring, suitable for O 167, O 173, O 188, O 190.



- Retro-fittable in upper baskets without magnetic spray arm

EAN: 4002514569857 / Material number: 07349850 / Old Material Number: 69113501D

Product classification	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Magnetic spray arm	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	75
External dimensions, net depth in mm	605
External dimensions, gross height in mm	47
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	615
Net weight in kg	0.31
Gross weight in kg	0.40