

MF 27/28-1

Transfer trolley

for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- Simple and convenient handling thanks to rounded handles
- Dockable or lockable to rolling conveyor etc. with both sides
- Infeed height 850 mm, -100 mm, +150 mm
- With removable drip tray H 70, W 603, D 866 mm

EAN: 4002514594026 / Material number: 07397640 / Old Material Number: 69200111D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Storage/transport system	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	1050
External dimensions, net width in mm	740
External dimensions, net depth in mm	930
External dimensions, gross height in mm	1010
External dimensions, gross width in mm	730
External dimensions, gross depth in mm	950
Net weight in kg	32.00
Gross weight in kg	45.00