

E 942/3

Module

for the optimum loading of up to 116 pipettes.



- Suitable for 116 pipettes up to 580 mm in height
- Compartment size ?? x ??

EAN: 4002514605036 / Material number: 07459390 / Old Material Number: 69594204D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert/module for pipettes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	116
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	279
External dimensions, net width in mm	558
External dimensions, net depth in mm	352
External dimensions, gross height in mm	420
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	570
Net weight in kg	35.00
Gross weight in kg	45.00