

E 498

Nozzle

for optimum loading of double and single tube breathing masks. Height 238 mm.

EAN: 4002514811673 / Material number: 07760120 / Old Material Number: 69749801D



- H 238 mm, with pipe 120 mm at the side

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	238
External dimensions, gross height in mm	80
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	250
Net weight in kg	0.12
Gross weight in kg	0.12