

E 900-5/2  
 Mobile unit  
 for the optimum cleaning of instruments on 5 levels.



- Loading capacity: 15 DIN mesh trays\*
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets ML)

\* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002514816418 / Material number: 07765760 / Old Material Number: 69590006D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	15
SPRI mesh trays [number]	10
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	605
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	680
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	650
Net weight in kg	12.85
Gross weight in kg	12.85
<b>Companion pack</b>	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1