

## DGG 11

Unperforated steam cooking container  
for cooking food in gravy, stock, water (e.g. rice, pasta).

EAN: 4002514177335 / Material number: 08234900 / Old Material Number: 22996174D



- Exterior dimensions W x D: 325 x 265 mm
- Exterior dimensions W x D x H: 325 x 354 x 20 mm
- Made from stainless steel
- Made from treated stainless steel, wall thickness 0.6 mm
- For steam oven container DGG 3
- Professional accessories manufactured in Europe
- Not suitable for Miele steam pressure ovens