

E 417/1

Insert

For the optimum loading of up to 30 ear and nose specula.

EAN: 4002514984605 / Material number: 09033260 / Old Material Number: 69541702D



- suitable for approx. 30 ear and nose specula
- 280 compartments approx. 10 x 10
- The mesh width on the base is 4 mm
- For use in upper and lower basket

|   |                 |
|---|-----------------|
| <b>Product category</b>                 |                 |
| Washer-disinfectors, medical            | •               |
| <b>Machine affiliation</b>              |                 |
| G 7831                                  | •               |
| PG 8535                                 | •               |
| PG 8536                                 | •               |
| PG 8562                                 | •               |
| PG 8582                                 | •               |
| PG 8582 CD                              | •               |
| PG 8592                                 | •               |
| <b>Product category</b>                 |                 |
| Insert for ENT instruments              | •               |
| <b>Product properties</b>               |                 |
| Material                                | Stainless steel |
| Colour                                  | Stainless Steel |
| Number of compartments                  | 280             |
| <b>Area of application</b>              |                 |
| Reprocessing of ENT instruments         | •               |
| <b>Standard electrical connection</b>   |                 |
| Number of phases                        | 0               |
| Voltage in V                            | 0               |
| Frequency in Hz                         | 0               |
| Total rated load in kW                  | 0.00            |
| Fuse rating in A                        | 0               |
| <b>Dimensions and weight</b>            |                 |
| External dimensions, net height in mm   | 80              |
| External dimensions, net width in mm    | 220             |
| External dimensions, net depth in mm    | 435             |
| External dimensions, gross height in mm | 85              |
| External dimensions, gross width in mm  | 225             |
| External dimensions, gross depth in mm  | 465             |
| Net weight in kg                        | 1.30            |
| Gross weight in kg                      | 1.72            |