

APH 591
Paper roll
for protocol printer APH 500 PRT100 and APH 510 PRT110.



- 5 paper rolls per pack
- Width: 76 mm mm

EAN: 4002515000038 / Material number: 09063410 / Old Material Number: 68817003D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Small sterilisers	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PS 1201 B	•
PS 1201 B EXCELLENCE	•
PS 1202 B	•
Product category	
Process documentation/operating documentation/communication	•
Product properties	
Colour	White
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of slides	•
Reprocessing of butyrometers	•
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•

APH 591

Paper roll

for protocol printer APH 500 PRT100 and APH 510 PRT110.



EAN: 4002515000038 / Material number: 09063410 / Old Material Number: 68817003D

Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of measuring cylinders	•
Reprocessing of containers, e.g. soda lime containers	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	79
External dimensions, net width in mm	76
External dimensions, net depth in mm	79
External dimensions, gross height in mm	105
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	210
Net weight in kg	1.39
Gross weight in kg	1.49