

ProCare Universal 61 - 6 kg  
 Dishwasher salt, 3 x 2 kg  
 For optimum regeneration of the integrated softener.



- Prevents calcium deposits and protects wash items

EAN: 4002515075319 / Material number: 09195780 / Old Material Number: 66916102EU1

<b>Product category</b>	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
<b>Applications</b>	
Reactivation of water softeners	•
<b>Dispensing recommendation</b>	
pH value at 20 °C	10,9
<b>Product properties</b>	
Medium type	Reactivation salt
Ingredients	Coarse special salt
Prevents limescale and protects load items	•
<b>Dimensions and weight</b>	
External dimensions, gross depth in mm	225.00
<b>Product properties</b>	
Medium type	Reactivation salt
Composition	Coarse special salt
Prevents limescale and protects load items	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	320
External dimensions, net width in mm	104
External dimensions, net depth in mm	226
External dimensions, gross height in mm	320
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	225
Weight in g	6000
Gross weight in kg	6
<b>Compliance, test certificates, and symbols</b>	
GHS/CLP compliant	•