


U 504
Basket
For the optimum loading of cutlery.



- With 9 plastic cutlery holders
- Compartment size 106 x 106 mm

EAN: 4002515250112 / Material number: 09508670 / Old Material Number: 67650401D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Basket for cutlery	•

U 504

Basket

For the optimum loading of cutlery.



EAN: 4002515250112 / Material number: 09508670 / Old Material Number: 67650401D

Product properties	
Material	Wire, rilsan-coated, Plastic
Colour	Grey, White
Number of compartments	9
Compartment width in mm	106
Compartment length in mm	106
Dimensions and weight	
External dimensions, net height in mm	145
External dimensions, net width in mm	500
External dimensions, net depth in mm	500
External dimensions, gross height in mm	260
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	530
Net weight in kg	2.24
Gross weight in kg	3.21